Product Description

Telio® Link is a dual-curing (light- and self-curing) temporary composite cement to be used for the esthetic temporary cementation of provisional crowns, bridges, inlays, onlays and veneers. Telio Link is compatible with temporary restorations fabricated with resin-based crown and bridge materials (e.g. Telio C&B) as well as with restorations made of Telio CAD.

Shades

Transparent, A3

Composition

Telio Link is composed of bismethacrylates (56 wt.%) and fillers (43 wt.%). Initiators, stabilizers and pigments are additional ingredients.

Indications

- Temporary cementation of provisional restorations.
- The product has been designed to remain in the mouth for a maximum of 6 weeks.

Contraindication

The product should not be used in patients with a known allergy to any of its ingredients.

Side effects

In rare cases, components of Telio Link may lead to a sensitizing reaction in susceptible individuals. In these cases, the use of the product should be discontinued.

Interactions

The eugenol-free formula makes Telio Link an ideal temporary cement even if adhesive luting techniques are subsequently employed. Eugenol-containing cements are contraindicated in these cases, as they may affect the polymerization of the luting composite used for the permanent cementation of the final restoration.

To prevent Telio Link from permanently adhering to composite-based core build-up material, freshly placed and prepared composite cores should be coated with a thin layer of separator. Water-soluble glycerine gel (e.g. Liquid Strip) may be used for this purpose. If dentin adhesives or protective varnishes are used after preparation, subsequent isolation with petroleum jelly or glycerine gel is required.

Application

To reduce dentinal hypersensitivity and postoperative sensitivity, Telio Desensitizer should be applied according to the Instructions for Use of Telio Desensitizer before temporary cementation.

Telio® Link**Telio® Link**

EN Instructions for Use	– Dual-curing resin based temporary cement	HR Upute za uporabu	– Dvostruko polimerizirajući privremeni cement na bazi smole	RO Instrucțiuni de utilizare	– Ciment temporar pe bază de rășină cu dublă polimerizare
DE Gebrauchsinformation	– Dualhärtender temporärer Befestigungszement auf Compositebasis	CS Návod k použití	– Duálně tuhnoucí provizorní cement na bázi pryskyřice	UK Інструкція щодо використання	– Тимчасовий цемент подвійної полімеризації на основі смоли
FR Mode d'emploi	– Colle provisoire duale	SK Návod na používanie	– Dočasný cement na báze duálne vytvrdzovanej živice	ET Kasutamisinstruktsioonid	– Ajutisel tsemendil põhinev kaksikkõvastuv vaik
IT Istruzioni d'uso	– Cemento composito provvisorio a polimerizzazione duale	HU Használati utasítás	– Kettős keményedésű gyantaalapú ideiglenes cement	LV Lietošanas instrukcija	– Divkārtšā cietināšanas sveķu bāzes pagaidu cements
ES Instrucciones de uso	– Cemento provisional en base a composites de polimerización dual	SR Упутство за употребу	– Привремени цемент на бази смоле за дуалну полимеризацију	LT Naudojimo instrukcija	– Dvigubo kietėjimo devos pagrindo laikinasis cementas
PT Instruções de Uso	– Cimento temporário, baseado em resina, com polimerização dual	MK Упатство за употреба	– Привремен цемент на база на двојно зацврстувачка смола		
FI Käyttöohjeet	– Kaksoskovetteinen yhdistelmämuovi-pohjainen väliaikaissementti	NO Bruksanvisning	– Dualhærdende provisorisk sement på komposittbasis		
NL Gebruiksaanwijzing	– Duaal uithardend, tijdelijk bevestigingscement op composietbasis	EL Οδηγίες Χρήσεως	– Ρηττωδής κοιλία διπλού πολυμερισμού για προσωρινή συγκόλληση		
		TR Kullanma Talimatı	– Çift sertleşen rezin tabanlı geçici siman		
		RU Инструкция по применению	– Временный фиксирующий цемент двойного отверждения на композитной основе		

Application of the Telio Desensitizer

1. The dentin surfaces should be dry and clean.
2. Apply Telio Desensitizer on the dentin and brush into the dentin for approximately 10 seconds using a suitable instrument (brush, application brush).
3. Carefully disperse excess with blown air / dry with an air syringe.

Application of the Telio Link

1. Remove the cap of the double-push syringe by turning it 1/4 turn counter clockwise (discard the cap, do not re-use it!) and replace it with a mixing tip. If the double-push syringe has been used before, replace the used mixing tip with a new one just before dispensing Telio Link.
2. When attaching the mixing tip to the syringe, make sure that the guidance of the double-push syringe is aligned with that of the mixing tip. Push the tip completely down until the notch on the mixing tip is aligned with that of the double-push syringe. Secure the mixing tip in place by gripping the coloured base and turning it 1/4 turn clockwise. The double-push syringe contains predosed amounts of Telio Link base and catalyst which are automatically mixed and dispensed when the two components are extruded. Thus the material may be directly applied to the provisional restoration.
3. Spread a layer of Telio Link over the dry inner surfaces of the provisional restoration and/or the prepared, cleaned and dried tooth surfaces. The working time is approx. 2.5 to 3 minutes at a room temperature of 23 °C / 73 °F. Incorporation of Telio CAD temporary restorations: Abrade the inside surfaces of the restoration by abrasive (with Al₂O₃, 100 µm grit at 1–2 bar/15–29 psi) or with a coarse diamond bur.
4. Seating of the provisional restoration and removal of excess cement: Seat the provisional restoration with firm but light pressure. Excess cement can be removed using different methods.
 - a) Removal of excess with additional light-activation (quarter technique) Light-cure excess material with the polymerization light (e.g. Bluephase® Style,) for 2–4 seconds per quarter surface (mesio-oral, disto-oral, mesio-buccal, disto-buccal) at a distance of approx. 0–10 mm. Following this, the excess material can be easily removed with a scaler as it features a viscoplastic consistency. Subsequently, light-cure all margins again for 10 seconds (>1000mW/cm²; e.g. Bluephase Style).
 - b) Removal of excess in the unpolymerized state Remove uncured excess material immediately after seating using a microbrush / brush / foam pellet / dental floss or scaler. Subsequently, wait until the self-curing process is finished (approx. 3 minutes) or accelerate the polymerization process by light-curing the cement layer with a curing light (>1000mW/cm²; e.g. Bluephase Style) for 10 seconds per aspect.
 - c) Removal of excess in the polymerized state Approx. 3 minutes after having seated the restoration, remove excess material carefully with a scaler or other instruments.

Leave the used mixing tip on the syringe.

It will serve as a seal until the next use!

Special notes

- To ensure an easy clean-up and complete removal of excess, we recommend holding back the lip or cheek surfaces during the setting process, to prevent the excess material from spreading out into a thin layer.
- During the self-curing process, the patient should not exert any pressure on the provisional restoration.

Working time

Approx. 2.5 to 3 minutes at a room temperature of 23 °C / 73 °F.

Self-curing time

Four minutes at a temperature of 37 °C / 98 °F.

Exposure to light will accelerate the curing process.

Note

- Extruding material without the mixing tip in place may result in a cross-contamination of the Telio Link components.
- Like all other composites, Telio Link is subject to oxygen inhibition. Therefore, the behaviour of the cement when mixed on a mixing pad outside the mouth cannot be taken as a reference for the intraoral setting behaviour.
- While the temporary restoration is in place, the material may undergo some change in colour. However, this does not affect the functionality.

Precautions

- Prevent any contact of uncured material with the skin/mucous membrane and eyes. Unpolymerized Telio Link may cause slight irritation and in very rare cases, sensitization may occur if the skin is exposed repeatedly.
- Commercially available medical gloves do not offer protection against the sensitizing effect of methacrylates.

Storage and shelf life

- Shelf life: see expiration date on label/packaging
- Do not use Telio Link after the indicated date of expiration.
- Storage temperature: 2 to 8 °C / 36 to 46 °F
- Only use Telio Link at room temperature (cooled material is more viscous and exhibits a slower setting behaviour).

Keep out of the reach of children!

For use in dentistry only!

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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