

An extraordinary family of curing lights.



Step up to a new class of curing lights  
with the Bluephase G2 and Bluephase 20i


**Bluephase**<sup>®</sup>  
LED CURING LIGHTS

**ivoclar**  
**vivadent**<sup>®</sup>  
passion vision innovation

# Bluephase® - an LED for every use

With its specially developed **Polywave® LED**, the 2<sup>nd</sup> generation of Bluephase curing lights sets new standards in the dental practice.

In contrast to conventional LED devices, new Polywave LED achieves an optimal broadband range from 385 to 515 nm, which is similar to the spectrum of halogen lights that served as the benchmark. The Polywave LED light is therefore suitable for all types of photoinitiators.



## Every material due to Polywave® LED

Due to the two different LEDs that are used – one with approx. 410 nm and the other approx. 470 nm dominant wavelength – Bluephase is suitable for all light-curing dental materials.



## Every indication due to continuous cooling

The internal fan allows for continuous operation without any clinical limitations – Bluephase allows for the placement of multiple indirect restorations without interruptions.



Simply remove power cord from the base.



Attach to handpiece and continue working.

## Every time due to Click & Cure


The proven **Click & Cure** function enables users to avoid irritating waiting times in the event the state-of-the-art lithium-polymer battery is depleted. The handpiece can be connected to the power cord of the charging base with just one click.

# Bluephase® G2

## extraordinary reliability

The battery-operated **Bluephase G2** featuring a light intensity of **1,200 mW/cm<sup>2</sup>**, three proven programs and a rotatable light probe with a constant diameter of 10 mm represents the new Bluephase market standard. Awards from **renowned institutions** are a testament to success and its perfected quality. The high light intensity and the design specific light collimation enables Bluephase to cure deep areas while offering the shortest curing times possible.

Ivoclar Vivadent provides the best prerequisites for a durable esthetic composite restoration by offering optimally coordinated products and outstanding precision in terms of light intensity. This is also supported by a field study conducted by the American Dental Association (ADA).



Source: R. Price, Dalhousie University Halifax, 2007  
 \* Not registered trademarks of Ivoclar Vivadent AG.

### Field test on light intensity (mW/cm<sup>2</sup>)

	Value indicated by manufacturer	Mean value measured
Bluephase**	1,100 (± 10 %)	1,076
L.E. Demetron II*	1,300	1,010
Coltolux LED*	1,400	609
SmartLite IQ*	900	544
FlashLite 1401*	1,400	439

Source: ADA, American Dental Association, Professional Product Review, Vol. 1 Issue 2, autumn 2006 \* Not registered trademarks of Ivoclar Vivadent AG.  
 \*\* Original Bluephase (not G2).

Light intensity of production curing lights was measured against manufacturer's specifications. Ivoclar Vivadent uses only high grade components and ensures a very narrow tolerance band through rigorous quality inspections of components and finished products.



red dot design award  
honourable mention 2008



# Bluephase® 20i

extraordinary power

The **battery-operated Bluephase 20i** combines the maximum light intensity of **2,000 mW/cm<sup>2</sup> in the Turbo program** with an extremely short curing time of **5 seconds** for light and dark composites while being gentle to the pulp and the soft tissue. Utilizing polywave broadband range output for curing of all common materials with confidence.

The full capacity of the Bluephase 20i is particularly useful when consistent and maximum performance is required, for instance, when indirect all-ceramic restorations are placed or orthodontic brackets are bonded. Due to the high efficiency, high output and the integrated fan for continuous cooling, adhesively cemented IPS Empress® and IPS e.max® restorations are polymerized in no time.



Four easy-to-use programs – Turbo for maximum performance, High Power for rapid curing, Low Power for curing areas near the pulp and Soft Start for stress-reduced polymerization.












The high light intensity allows polymerization in the shortest possible time in every clinical situation.



**Short curing times of just 5 seconds**  
With broadband range for curing of all common materials with confidence.

# Bluephase® – Technical data



	Bluephase® G2 1,200 mW/cm <sup>2</sup> +/-10%	Bluephase® 20i 2,000 - 2,200 mW/cm <sup>2</sup> LED Class 2	Bluephase® meter 300 - 2,500 mW/cm <sup>2</sup> +/- 20%
<b>Every material</b> (wavelength range)	✓ 385 - 515 nm	✓ 385 - 515 nm	✓ 385 - 515 nm
<b>Every indication</b> (continuous operation for at least 10 min)	✓	✓	Measuring the light intensity of LED curing lights
<b>Every time</b> Click&Cure (optional corded operation)	✓	✓	
<b>Maximum curing time for composites</b>	15 sec	10 sec	
<b>Maximum curing time for Ivoclar Vivadent composites</b>	10 sec	5 sec	
<b>Curing programs</b>			
<b>TURBO</b>  	—	2,000 mW/cm <sup>2</sup>	
<b>HIGH Power</b>  	1,200 mW/cm <sup>2</sup>	1,200 mW/cm <sup>2</sup>	
<b>LOW Power</b>  	650 mW/cm <sup>2</sup>	650 mW/cm <sup>2</sup>	
<b>SOFT Start</b>  	650 / 1,200 mW/cm <sup>2</sup>	650 / 1,200 mW/cm <sup>2</sup>	
<b>Light probe</b>	10 mm black	10>8 mm black	
<b>Power supply</b>	Lithium-polymer battery capacity: approx. 60 min charging time: approx. 2h	Lithium-polymer battery capacity: approx. 45 min charging time: approx. 2h	3 x LR6 AA 1,5 VDC
<b>Display</b>	OLED color display	OLED color display	Digital LCD display
<b>Warranty</b>	3 years (battery 1 year)	3 years (battery 1 year)	3 years

Technical data

# Bluephase® – Delivery forms and accessories



## 607920 Bluephase® G2 kit

- 1 - handpiece
- 1 - 10mm light probe
- 1 - battery
- 1 - charging base
- 1 - power supply
- 50 - protective sleeves
- 3 - Anti-glare cones

## 607921 Bluephase® G2 kit with power meter

- 1 - Bluephase® G2 kit
- 1 - Bluephase power meter

## Bluephase® G2 accessories

- 607922 Bluephase power meter
- 608532 Handpiece (with battery and 10mm light probe)
- 608535 Battery
- 608537 Light probe 10mm
- 608538 Light probe 6/2 mm Pin-Point
- 608554 Protective sleeves (5 x 50 pieces)
- 551756 Anti-glare cones (3 pieces)
- 592496 Anti-glare shield



## 613735 Bluephase® 20i kit

- 1 - handpiece
- 1 - 10/8 mm light probe
- 1 - battery
- 1 - charging base
- 1 - power supply
- 50 - protective sleeves
- 1 - Anti-glare shield
- 3 - Anti-glare cones

## 613752 Bluephase® 20i kit with power meter

- 1 - Bluephase® 20i kit
- 1 - Bluephase power meter

## Bluephase® 20i accessories

- 607922 Bluephase power meter
- 613753 Handpiece (with battery and 10/8mm light probe)
- 627300 Battery
- 627389 Light probe 10/8mm
- 608538 Light probe 6/2 mm Pin-Point
- 608554 Protective sleeves (5 x 50 pieces)
- 551756 Anti-glare cones (3 pieces)
- 592496 Anti-glare shield

*Every material. Every indication. Every time.*

100% CUSTOMER SATISFACTION  
GUARANTEED!

[ivoclarvivadent.com](http://ivoclarvivadent.com)

©2010 Ivoclar Vivadent, Inc. Bluephase is a registered trademark of Ivoclar Vivadent, Inc. Flashlight 1401, Coltolux LED, Demi, L.E. Demetron II, Smartlite IQ2, Smartlite PS and Elipar Freelight 2 are not registered trademarks of Ivoclar Vivadent, Inc.  
627880 Rev. 3 12/2012

### United States

175 Pineview Drive  
Amherst, NY 14228  
800.533.6825  
Fax: 716.691.2285  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

### Canada

1 - 6600 Dixie Road  
Mississauga, Ontario L5T 2Y2  
800.263.8182  
Fax: 905.670.3102  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**ivoclar**  
**vivadent**  
passion vision innovation