# Email Preparator blue / blau

## Instructions for Use

### Gebrauchsinformation

### Mode d'emploi

### Istruzioni d'uso

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### Bruksanvisning

### Brugsanvisning

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### Productinformatie

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For dental use only!

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

CE 0123

558172/0805/WE3/H

English

## Instructions for Use

### Description

Email Preparator blue is a gel for enamel etching and dentin conditioning. It is used in conjunction with sealings and direct composite restorations, as well as lab-fabricated all-ceramic or composite restorations placed with the adhesive technique .

### Composition

Email Preparator blue contains phosphoric acid (37 wt.% in water), polyvinyl alcohol, and pigments.

### Indication

Email Preparator blue is used in the enamel etch or total etch techniques for the following indications:

- Composite restorations
- Adhesive cementation of inlays (composite and all-ceramic), crowns, veneers, and adhesive bridges
- Splints
- Brackets
- Fissure sealings

### Contraindication

If the stipulated working technique cannot be used, the material should not be applied.

### Side effects

None known to date.

### Interaction

Glass ionomer-based bases/liners may be dissolved by Email Preparator blue.

### Application procedure

Prior to applying the material, fix the enclosed cannula to the top of the bottle and bend it according to requirements. Apply Email Preparator blue on the enamel and then on dentin and allow a reaction time of 15 seconds. After that, thoroughly rinse with water spray and dry the tooth surfaces with oil-free air. Avoid dehydrating the dentin. The etched enamel surface should have a chalky white appearance. If this is not the case, repeat etching. Should

ivoclar  
vivadent<sup>®</sup>  
clinical

the tooth surfaces accidentally be contaminated (e.g. with saliva), the etching procedure has to be repeated.

The reaction time on unprepared enamel (e.g. with fissure sealings) is 30–60 seconds. If the enamel is selectively etched, Total Etch has to be left to react for 15–30 seconds.

#### Warning

Email Preparator blue contains phosphoric acid and is corrosive. Avoid contact with the eyes, mucous membrane and skin (goggles are recommended for both patients and practitioners).

After contact with the eyes, immediately rinse with copious amounts of water and see a physician/ophthalmologist. After skin contact, immediately rinse with soap and water.

#### Storage

- Do not use Email Preparator blue after the indicated date of expiration.
- Store at 2–28 °C
- Date of expiration: see note on bottle/packaging.
- By repeated cooling of the etching gel to refrigerator temperature or below 0 °C / 32 °F, a cryogel may develop. Thus, the viscosity increases. This process, however, is reversible (e.g. by heating up the gel to 50–80 °C / 122–176 °F for 5–10 minutes) and does not damage the product.

**Keep out of the reach of children.**

**For use in dentistry only.**

**Date information prepared:** 03/2004

#### Manufacturer

Ivoclar Vivadent AG

FL-9494 Schaan / Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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