EN INSTRUCTIONS FOR USE

MODELLATION

Design the framework in a reduced anatomic shape taking the planned veneer into consideration... Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm.

SPRUNG

Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center.

INVESTMENT

Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy).

BURN-OUT

The suggested burnout temperature: 800-850C/1470-1560F

MELTING AND CASTING

Use a separate ceramic crucible for each alloy and preheat the crucible in the burnout furnace. Always use fresh alloy. Depending on the type of casting machine, follow the manufacturer's instructions for use.

Casting Temperature: 1240-1350C/2265-2460F

METAL PREPARATION

Carefully divest and clean the object with Al₂O₃. Do not use a hammer for divesting the object to prevent deformation.

OXIDATION

Blast surface with 50-100 micron Al₂O₃ at max. 5.2 bar/75 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam.

SOLDERS AND FLUXES

Design the soldering pattern(s) as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip.

Pre-Solder: Special High Fusing White Ceramic Solder, High Fusing White Ceramic Solder, High Fusing Bondal Flux

Laser weld material: Laser Ceramic White

POLISHING

After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.

INDICATIONS

Recommended for short span bridges, telescope crowns, conus crowns, posts, long span bridges, implant superstructures, partial dentures, PFM crowns.

CONTRAINDICATIONS

For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS

In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS

Galvanic effects may occur between different alloys in the same oral environment.

For additional information look into the alloy property chart.

CONTRAINDICATIONS

For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS

In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS

Galvanic effects may occur between different alloys in the same oral environment.

For additional information look into the alloy property chart.

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