ivoclar
vivident[®]
technical

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Litho USA 5/05
Rev. 8/06

Caution: US Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in U.S.A. unless otherwise indicated

Made in Canada (fabrique au Canada)

LOT

www.ivoclarvivadent.com

Quality Management System Certified

ivoclar vivident technical

ivoclar
vivident[®]
technical

Recommended for onlays, 3/4 crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

Au	Zn	Ga	Ru	Ir	Re	Other
2.5	1.5	6.0	<1.0	<1.0	<1.0	<1.0

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Litho USA 5/05
Rev. 8/06

Caution: US Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in U.S.A. unless otherwise indicated

Made in Canada (fabrique au Canada)

LOT

www.ivoclarvivadent.com

Quality Management System Certified

ivoclar vivident technical

INSTRUCTIONS FOR USE

EN

MODELLATION
Design the framework in a reduced anatomic shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.

SPRUNG
Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

INVESTMENT
Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instructions.

BURN-OUT
The suggested burnout temperature: 750-820C/1380-1510F

MELTING AND CASTING
Use a separate ceramic crucible for each alloy and preheat the crucible in the burnout furnace. Used and new alloys must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturer's instructions for use. If you use Williams Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting bench cool to room temperature.

Casting Temperature: 1390-1450C/2535-2645F

METAL PREPARATION
Carefully divest and clean the object with AL₂O₃. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments.

OXIDATION
Blast surface with 50-100 micron AL₂O₃ at max. 5.2 bar/75 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam. Place the object on the firing tray and provide adequate support. Place the tray in a porcelain furnace set at a low temperature of 650C/1200F and increase the furnace temperature to 1035C/1900F without vacuum with 5 min. hold time at top temperature.

HEAT TREATMENT
Hardening: 650C/1200F for 30 minutes; air cool.

SOLDERS AND FLUXES
Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.

Pre-Solder: High Fusing White Ceramic Solder, High Fusing Bondal Flux
Post-Solder: .650 Fine Solder, .615 Fine Solder, Low Fusing White Gold Solder, Bondal Flux
Laser weld material: Laser Ceramic White

POLISHING
After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.

INDICATIONS
Recommended for onlays, 3/4 crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

CONTRAINDICATIONS
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS
In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS
Galvanic effects may occur between different alloys in the same oral environment.
For additional information look into the alloy property chart.

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