span bridges, PFM crowns, implants, 3/4 crowns, telescope crowns, conus crowns, posts, short and long
recommendation for onlays, 3/4 crowns, telescope crowns, implants superstructures, partial dentures

Instructions
Composition
Alu. Pu. Pd. Ag. Cu. S. Zn. In. Ga. Ru. Ir. Re. Other
- 37.0. - 55.9. 1.0. - - - - -

Kayttoohjeet
Olytöt ja yhdysk
Bruksanvisnings
Instruktioner de
bruksanvisning
Bruksanvisnings
Produktinformations

Austria
Ivoclar Vivadent GmbH
Liechtenstein
Ivoclar Vivadent AG
Benderstrasse 2
FL-9494 Schaan
Tel. +423 235 35
Fax +423 235 30
www.ivoclarvivadent.com

Australia
Ivoclar Vivadent Pty. Ltd.
5/5 Overses Drive
P.O. Box 367
Noble Park Vic 3174
Tel. +61 3 979 595 99
Fax +61 3 979 595 46
www.ivoclarvivadent.com.au

ISO 13485
Quality Management System Certified
Made in U.S.A. unless
otherwise indicated
Made in Canada
(fabrique au Canada)

USA
Ivoclar Vivadent, Inc.
175 Pineview Drive
Amherst, NY 14228
Tel. +1 800 522 6826
Fax +1 716 691 2285
www.ivoclarvivadent.us.com

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White Ceramic Alloy
SI@93

SPRUING METHOD

RECOMMENDATIONS
DIRECT:
single crowns,
inlays and onlays

INSTRUCTIONS:

1. Select a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration.
2. Maintain the reservoir(s) within the heat center of the investment; keep the restoration(s) approximately 5 mm from the end of the investment and no closer than 5 mm from the sides.
3. Make sure that the sprue is connected to the thinnest part of the restoration.
4. The connection between the sprue and the restoration should be flared ("trumpet" configuration); eliminating turbulence of the alloy (causing investment erosion) and facilitating an undisturbed flow of alloy during casting and solidification.
5. Maintain proper alloy weight in order to eliminate the negative effect of a too large button during solidification. The rule of thumb for proper alloy weight determination is: total wax weight x specific gravity of the alloy = proper alloy weight.

SUGGESTIONS:
1. Thickness and design of wax pattern: follow your manufacturers instructions.
2. Use chill-vents when casting heavy and/or large restorations.

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INDIRETTO:

multiple units
e diverse
single crowns

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