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Instructions for use

1 Waxing up

Minimum wall thickness of 0.3 mm for single crowns and minimum of 0.4 mm for bridge abutment crowns. Cut back frameworks anatomically to allow for the planned veneer and to ensure stable shaping of frameworks for larger span widths; avoid sharp transitions. Create large-area surfaces for the planned soldering.

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9 Aushärtung

(when used as a casting alloy or with composite/resin veneers)

Observe the instructions of the veneering material manufacturers.

2 Spruing

For single crowns and bridges of up to max. 3 units, direct spruing (min. 3 mm diameter) or the use of pear-shaped sprues is recommended. In general, the cast bar (approx. 4–5 mm diameter) is fitted with guide channels (3–4 mm diameter) and connecting channels (approx. 3 mm x approx. 3 mm diameter) and if necessary cooling ribs (1 mm diameter). The wax pattern including the sprues must be weighed in order to determine the required amount of alloy (see wax conversion table: wax weight x density = amount of alloy in g).

3 Investing

Line the casting flask with a suitable ceramic liner. Mix the recommended investment material (see data table (A)) under vacuum and invest (follow the manufacturer's instructions carefully).

4 Preheating

With conventional gradual preheating methods, hold the first preheating step at 280 °C / 536 °F for 30–60 min; further preheating according to the instructions of the investment material manufacturer. Once the final temperature is reached (see data table (B)), the holding time, depending on the size of the flask, is as follows: 20 min (1x), 40 min (3x), 50 min (6x), 60 min (9x).

When using speed preheating methods, it is essential to follow the instructions of the investment material manufacturer.

5 Melting, casting

See data table (C, D, E) for casting crucible, casting temperature and post-melting time, depending on the method used. It is advisable to use a separate crucible for each alloy. Use only very clean casting residues. The recommended ratio of used material to new material is 1:1.

6 Finishing, cleaning

After slow cooling until warm to the touch, divest and if appropriate blast with glass beads (approx. 50 µm) or dip (Wilavest Sonic Liquid). Finish framework with sharp carbide cutters or ceramic-bonded grinding instruments at minimal pressure. Then blast surfaces with pure aluminium oxide (approx. 110 µm) at 2 bar pressure. Subsequently, steam-clean the framework thoroughly and if necessary degrease with ethyl alcohol.

7 Oxidizing

■ Side effects: In rare cases, sensitivity or an allergic reaction to one or several components of the alloy may occur.

■ Contraindications: If a patient is known to be allergic or sensitive to one or several of the alloy ingredients, a physician should be consulted.

■ Interactions: Several different types of alloys in the oral cavity may cause galvanic reactions.

8 Soldering, laser welding

For soldering before and/or after ceramic firing and before resin veneering, use solders as stipulated in the data table (H); soldering gap of 0.05 to 0.2 mm.

Allow the soldered object to cool slowly after soldering.

For laser welding, use welding wire as stipulated in the data table (J).

Caution: Pay special attention to the weld seam geometry!

9 Hardening

(when used as a casting alloy or with composite/resin veneers)

The alloy hardens by itself during the slow cooling after casting or after soldering (s.). If necessary, maximum hardening (v) can be achieved by additional annealing (temperature/time, see data table (G)).

10 Veneering

All commercially available low-fusing, high-expanding ceramic materials with a CTE between 16.1 and 16.7 x 10⁶ K¹ and with a matching firing temperature (generally, min. 100 °C / 212 °F below the solidus temperature of the alloy) are suitable for ceramic veneering. A matrix firing in two stages is recommended; follow the instructions of the particular ceramic manufacturer for further ceramic processing; ensure adequate support.

Before composite veneering, dipping (Wilacid) may have to be carried out to remove oxide and flux residue after soldering and/or hardening processes.

Additional finishing or cleaning should be performed as described in point 6.

It is essential to observe the handling recommendations of the veneering material manufacturers.

11 Dipping, polishing

Remove residual oxide and any flux residues by dipping (Wilacid) or by blast polishing. Smooth the metal surfaces finished with carbide cutters using ceramic-bonded abrasives. Then polish to a semi-gloss shine with rubber polishers. Polish to a high-gloss finish with suitable polishing paste (Polirol U). Subsequently, clean using ultrasonic cleaning equipment or by careful steam blasting.

12 Preheating

With conventional gradual preheating methods, hold the first preheating step at 280 °C / 536 °F for 30–60 min; further preheating according to the instructions of the investment material manufacturer. Once the final temperature is reached (see data table (B)), the holding time, depending on the size of the flask, is as follows: 20 min (1x), 40 min (3x), 50 min (6x), 60 min (9x).

When using speed preheating methods, it is essential to follow the instructions of the investment material manufacturer.

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