SR Ivocap® System

EN Instructions for Use

Page 2

DE Gebrauchsinformation

Seite 8

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Page 14

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Pagina 20

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Pagina 26

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Página 32

CE 0123

EN ISO 20795-1: 2013

Rx ONLY





Date information prepared:
598366/2016-11-28/WE1/Rev.2






Manufacturer
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclarvivadent.com




The logo for Ivoclar Vivadent, featuring a series of dots above the word 'ivoclar' and a series of dots above the word 'vivadent'. Below 'vivadent' is the word 'technical' in a smaller, italicized font.
ivoclar
vivadent
technical

List of parts / Teileverzeichnis
Liste des composants / Elenco particolari /
Despiece / Lista de partes

Step-by-step procedure / Verarbeitung / Mise en oeuvre / Istruzioni d'uso / Instrucciones de Uso / Procedimento

Please read these instructions carefully to familiarise yourself with the operation.

Product description

The SR Ivocap System is an injection technique, which has been especially developed to compensate curing shrinkage. Heat/pressure curing allows shrinkage of the acrylate during polymerization to be compensated by inflowing material. The SR Ivocap injection technique, therefore, permits the fabrication of high quality products with excellent physical properties.

Composition

SR Ivocap High Impact

Powder:

Polymethyl methacrylate, copolymer, benzoyl peroxide, pigments

Liquid:

Methyl methacrylate, dimethacrylate (cross-linking agent), copolymer

SR Ivocap Clear

Powder:

Polymethyl methacrylate, dibenzoyl peroxide

Liquid:

Methyl methacrylate, dimethacrylate (cross-linking agent), stabilizer

List of parts

SR Ivocap System

A) Pressure apparatus "ID2"

- A1 = Manometer
- A2 = Extrusion piston
- A3 = Vise shank
- A4 = Safety loop

B) Clamping frame

- B1 = Ratchet
- B2 = Pressure ring with lever
- B3 = Pressing plate

C) Flask

- C1 = Flask half M
- C2 = Flask half W
- C3 = Flask cover
- C4 = Injection funnel

D) Capsule plunger

E) Investment aid

F) Funnel

G) Timer

Indication

SR Ivocap High Impact, SR Ivocap Clear

- Complete dentures
- Partial dentures
- Bases and relines
- Orthodontic appliances
- Bite guard splints

Contraindication

- Direct intraoral contact of unpolymerized material.
- If the patient is known to be allergic to any of the components listed.

Important notes

- Do not wet teeth with SR Ivocap monomer, as this will weaken the bond to the tooth
- Soak models thoroughly before the investment
- Thoroughly isolate inner flask surfaces before the investment
- Injection channels must have a diameter of 3–5 mm
- Thoroughly clean and isolate plaster surfaces
- Observe polymerization and cooling times precisely
- Contact of solvents or monomer with the polymerized denture base material can cause white discolouration
- If the cervical area is to be isolated, Ivoclar Vivadent recommends using A silicones with a Shore A hardness of >65 (e.g. Flexistone Plus/Virtual[®] Heavy Body Fast). Condensation-reaction silicones (C silicones) may have a negative effect on the denture base material. Other products may cause negative reactions of the denture base material.
- The monomer composition of SR Ivocap HI, SR Ivocap Clear and SR Ivocap Elastomer is not identical. Avoid a mix up of these three monomers as it causes product defects.

Warnings

- The monomer contains methyl methacrylate (MMA)
- MMA is highly flammable and irritating (flash point: +10 °C / 50 °F).
- MMA irritates the eyes, respiratory organs and skin
- Contact may evoke sensitive skin reaction.
- Avoid contact of skin with monomer and unpolymerized material. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Do not inhale vapours
- Keep the material away from sources of heat – do not smoke
- Do not dispose of the monomer in the drainage system
- Avoid electrostatic charging

Step-by-step procedure

① Investment (lower flask half)

Isolate the flask against plaster using Vaseline. Place the investment aid (1) in the flask, mix the plaster and invest the thoroughly water-soaked model. Remove excess plaster from the thermal insulating component (2) and the flask edges.

② Placing the funnel

Remove the investment aid after the plaster has set and position the injection funnel (3) and the funnel (4) in its place.

③ Injection channels

For complete upper dentures one injection channel is sufficient, fanning out from the funnel tip towards the palate (Fig 3a). In the case of complete lower dentures two channels to the dorsal extremities of the dentures are required. (Fig 3b). For cast metal frameworks with individual extension bases, either each base must be provided with a channel at the distal end or an interconnecting ring channel (dotted line) can be fabricated (Fig 3c). Prior to investing the antagonist part, the injection channels must be formed with pink wax. This can also be done after boil-out or before isolating by simply milling them in.

IMPORTANT:

- Injection channels must have a diameter of 3–5 mm
- Choose the shortest path
- Make sure not to damage the investment edges along the injection funnel; they form a seal towards the posterior shoulder of the flask.

④ Investment (upper flask half)

Place the upper flask half on the lower one after isolating it with Vaseline or Ivoclar Vivadent Separating Fluid. Mix the plaster and pour over teeth, avoiding the formation of bubbles, until the incisal edges and occlusal areas are lightly covered.

To make deflasking easier, apply a paper insert as a separation layer (the insert should be moistened with water).

Then, fill the flask to the brim, put on the cover and completely press down by hand.

Do not close the flask with a press!

⑤ Boil-out procedure

Preheat the flask in very hot water for at least 5 minutes. Remove wax and thoroughly boil out with hot and clean water. Allow the flask halves to cool to room temperature. Remove adhering plaster residue from the flask edges and thermal insulation component.

⑥ Isolation

Isolate the still moist plaster surfaces twice with Ivoclar Vivadent Separating Fluid. Separating with a suitable spray gun is also possible. Avoid pooling, eg in the gingivobuccal fold or interdental areas. Do **not** rinse isolated areas with water.

⑦ Capsule preparation

Remove monomer container from the capsule and break open by twisting off the sealing cap.

⑧ Pour the monomer into the capsule. Re-insert the monomer container into the capsule base prior to mixing. Briefly shake the mixture.

⑨ Mount the capsule in the Cap Vibrator, secure it with a rubber band and mix for 5 minutes. The SR Ivocap denture base material should form a ball in the process. If this is not the case after 5 minutes of mixing, the material can be mixed for an additional 30 seconds. The failure of the material to form a ball does not influence the quality of the denture base material.

⑩ Remove the empty monomer container, eg with compressed air, place the capsule on the capsule plunger and press the contents upward with light rocking movements. Allow air to escape through the capsule opening.

IMPORTANT:

- Start mixing immediately following the addition of the monomer (5 min in the Cap Vibrator).
- Do not mix longer than 7 minutes, as the resin may become too warm, which increases the risk of premature polymerization.
- There should be no air between the plunger and the contents.
- Previously mixed capsules can be refrigerated in the dark for 5 days. Allow them to regain room temperature before reuse.

⑪ Clamping of the flask

Place the injection funnel with the inserted funnel in the lower flask half. Carefully bring together the two flask halves and ensure the correct fit of the flask covers. Then, completely insert the flask into the clamping frame. Apply 3 tons / 6000 lbs pressure to the clamping frame with the flask in a hydraulic press (corresponds to about 80 bar / 1133 psi hydraulic pressure). At the same time, push the ratchet (B1) on the clamping frame to the right. After that, cease applying pressure and remove the clamping frame with the flask from the press.

⑫ Inserting the SR Ivocap capsules

Remove the cover from the capsule and fully insert the capsule into the flask.

⑬ Mounting the pressure apparatus

Fully extend the plunger of the pressure apparatus, e.g. by pulling the piston rod when the compressed air locking valve is open.

- ⑭ Place the SR Ivocap pressure apparatus on the flask, lift the safety loop (A4) and press the plunger into the capsule until the vise shanks are correctly positioned on the flask. Then pull down the safety loop until it snaps into place. Make sure that the plunger does not get jammed in the capsule.

⑮ Injection

Connect the pressure apparatus **with the compressed air locking valve closed** to the compressed air supply (6 bar / 85 psi). Slowly open the locking valve. The plunger will descend and press the SR Ivocap material into the mould.

IMPORTANT:

- **Fully extend the plunger of the pressure apparatus.**
- **Leave the material under pressure for 5 minutes before placing it in a water bath.**
- **SR Ivocap material which is left over or which has been in a mixed state for some time requires a longer injection time of 10 minutes.**
- **Allow leftover material to regain room temperature before use.**
- **When the red ring on the piston rod becomes visible during injection, mix a new capsule and replace the empty one.**

⑯ Polymerization

Place the SR Ivocap assembly in a suitable polymerization bath. The surface of the water should be covered with plastic floaters in order to prevent unnecessary loss of heat (isolation). When immersing the SR Ivocap assembly in the polymerization bath, make sure that no plastic floaters are trapped beneath the clamping frame.

- ⑰ The water level must reach the red mark on the clamping frame and must be maintained at this level throughout the polymerization period.

- The temperature of the water bath must be set in such a way that the water boils during the entire polymerization period.
- The polymerization time once the water **starts to boil** is exactly 35 minutes.
- Do not interrupt the boiling procedure by introducing other flasks.
- Content of residual monomer according to EN ISO 20795-1 < 2.2%. In order to reduce the content of residual monomer below 1%, the material must be polymerized for 90 minutes in boiling water. The residual material will set in the capsule during this time.

Consequences of incorrect water level during polymerization (Fig. 17)

over the red mark

- Capsule reservoir polymerizes; constant flow of the resin during polymerization becomes impossible.
- Waste of material due to polymerization of leftover material in the capsule.

below the red mark

- Partially unpolymerized areas in palatal, post dam or retromolar areas

⑱ Cooling

After completion of the polymerization procedure, remove the SR Ivocap assembly from the boiling water and immediately cool in **cold** water.

During the first 20 minutes of the cooling phase, the pressure in the clamping frame and the injection apparatus must remain unchanged.

After 20 minutes, the pressure apparatus can be removed. The clamping frame together with the flask, however, must remain in cold water for an additional 10 minutes.

Total cooling time: 30 minutes

- ⑲ A locking valve on the pressure apparatus maintains the injection pressure at a constant level during the cooling phase, irrespective of the compressed air supply. This results in increased mobility, since the cooling bath does not necessarily have to be near the air supply.

⑳ Deflasking

Mount the clamping frame in the press and re-apply 3 tons / 6000 lbs pressure. Move the ratchet (B1) back to the left and release pressure. Remove the flask from the clamping frame.

㉑ Recommended procedure

- Remove both plastic covers, insert flask with the narrow side to the front into the press. Place any suitable object between the plaster investment and the pressing plate, not touching the flask edges. Close the press with minimal pressure.
- Introduce a broad screwdriver or plaster knife into the joint beside the guide piece and lever the upper flask half upwards.
- Remove the flask from the press and repeat the same procedure with the flask reversed. The investment can now be removed from the flask without difficulty and without using force (e.g. hammer).
- Remove plaster in the usual way.
- Finish and polish as usual.

Repairs

All the SR Ivocap shades are also available in cold-curing form. Ivoclar Vivadent recommends its cold-curing polymers for repairs. The Instructions for Use of these materials provide detailed information about the repair methods available. Please see the individual Operating Instructions of the different apparatus.

Lagerungshinweise


- Store the material in a cool, dark and well-ventilated place
Unmixed: storage temperature of 2–28 °C / 36–82 °F
Mixed: max 5 days, cool and protected from light at 8 °C / 46 °F
- If the material is refrigerated, the capsule should be removed from the refrigerator one hour before use and brought to room temperature.
- Do not use the product after the indicated expiry date.
- Store out of the reach of children.

Classification according to EN ISO 20795-1


Type 1, Class 1

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.


Operating diagram

Wall attachment for
4 pressure apparatus
(on request)

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SR Ivocap High Impact

Pulver:

Polymethylmethacrylat, Co-Polymer, Benzoylperoxid,

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.