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Instructions for Use

Product description

Transparent, medium viscosity (Type 2), two-component addition silicone. The product sets at room temperature. Especially suitable for the fabrication of matrices (C&B technique) in dental laboratories.

Composition

Transil contains vinyl polysiloxane, polymethylhydrogen siloxane, organic platinum complex and silicate.

Applications

- Transparent matrix for the fabrication of Vectris frameworks
- Matrix in the C&B technique in conjunction with the wax up

Limitations of use

Transil should not be used intraorally.

Side effects

Not known.

Interactions

Latex gloves inhibit the curing of vinyl polysiloxanes. Therefore, Transil should not come into direct or indirect contact with latex gloves. For the processing of Transil, use either vinyl gloves or wash your hands thoroughly in advance.

Processing

Twist off the cap of the cartridge and mount the mixing tip. Subsequently, insert the cartridge into the dispenser. The base and catalyst are auto-mixed in a 1:1 ratio in the mixing tip. Extrude the material evenly until the matrix has been fabricated. After application leave the mixing tip on the cartridge, functioning as a cap. Use a new mixing tip for

each application. At room temperature (23 °C / 73 ° F), Transil features a processing time of 1–1.5 minutes and a setting time of 15 minutes.

Warning

After accidental contact with the eyes, immediately rinse with copious amounts of water. If the irritation continues, seek medical advice. If the material comes into contact with the skin, rinse with water and soap.

Storage

- Do not use Transil after expiry date.
- Store at 12-28 °C / 54-82 °F
- Expiry date: See label on packaging
- Keep away from sources of heat

**Keep out of the reach of children.
For use in dentistry only.**

Date information prepared

11/2003

Ivoclar Vivadent AG
FL-9494 Schaan / Liechtenstein

This material has been developed solely for use in dentistry. Processing should be done strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent S.A.
Ivoclar Vivadent AB
Ivoclar Vivadent, Inc.