
# Helioseal® F Plus

The efficient fissure sealant





Tight sealing of fissures with Helioseal® F Plus (Picture: M. Hauner, Ivoclar Vivadent AG,

# Helioseal® F Plus

Solve Spanish Spanish





Permanent molar before and after the application of Helioseal® F Plus.

# Precise dispensing

Helioseal® F Plus is a lightcuring, white-pigmented, fluoride-releasing fissure sealant. It comes in a syringe, which allows the product to be precisely dispensed and applied. With the fine tip the sealant is applied quickly and accurately without overfilling the fissure.

#### **Excellent flow properties**

Helioseal F Plus is flowable and it is easy to distribute on intricate surfaces – even in the upper jaw. The tight marginal seal provides protection against cariogenic germs. Fine filler particles ensure a smoothly sealed surface.

#### Fast curing

With the help of an appropriate curing light, such as Bluephase Style or Bluephase G4, Helioseal F Plus can be cured in only ten seconds. Due to its short light-curing time, precise application and excellent flow and wetting properties, the sealant offers fast and efficient treatment in children and adolescents.

#### Helioseal® F Plus - The efficient fissure sealant

## **Advantages**

New formulation featuring optimized consistency

Precise dispensing and application

Reduced exposure time

### **Benefits**

Optimum distribution on the intricate surfaces of pits and fissures

Efficient treatment

Fast, comfortable treatment for children



#### Helioseal® F Plus deliver forms



**Helioseal F Plus**Syringe containing
1.25 g

**Helioseal F Plus**Cavifil containing
0.1 g

Literature

Fischer K: Scientific Documentation Helioseal® F Plus, Ivoclar Vivadent



Production and distribution Ivoclar Vivadent AG Bendererstr. 2 9494 Schaan Liechtenstein Tel. +423 235 35 35 Fax +423 235 33 60 www.ivoclarvivadent.com Descriptions and data constitute no warranty of attributes. © Ivoclar Vivadent AG, Schaan/Liechtenstein. 709048/05.06.2019

