

SR Phonares® II Typ

SR Phonares® II Typ teeth are made of the same filled hybrid composite of the latest generation as the anterior teeth. The semi-anatomical mould is designed for the conventional setup method.



Advantages:

- High wear resistance
- Reduced basal design renders time-consuming reductions unnecessary
- Ideally suitable for covering implant abutments and other construction elements
- For the conventional setup method

Delivery forms

SR Phonares® II Typ Posterior Teeth

16 A–D, 4 Bleach shades
3 upper, 3 lower

8 teeth per set

639352

SR Phonares® II Lingual

SR Phonares® II Lingual teeth are made of the same filled hybrid composite of the latest generation that is used for the anterior teeth. The mould is designed for the lingualized setup method, which is why these teeth are the best choice for use in implantology.



Advantages:

- High wear resistance
- Reduced basal design renders time-consuming reductions unnecessary
- Ideally suitable for covering implant abutments and other construction elements
- For the lingualized setup method

Delivery forms

SR Phonares® II Lingual Posterior Teeth

16 A–D, 4 Bleach shades
3 upper, 3 lower

8 teeth per set

639351