

SR Orthotyp® S PE

The original posterior moulds are based on a detailed functional analysis of the stomatognathic system. The modern appearance meets high prosthetic and esthetic demands. In line with the principle of group function, the Orthotyp posterior moulds are set up in a one-to-two-tooth relation.



Advantages:

- Highly cross-linked DCL polymer
- Harmonious progression of crown lengths to create an esthetically pleasing transition from anterior to posterior moulds
- Wide tooth neck to hide esthetically placed structural elements
- Naturally designed surface structure to match the esthetic appearance of the anterior moulds

Delivery forms

SR Orthotyp® S PE Posterior Teeth

20 PE shades		
4 upper, 4 lower	8 teeth per set	657193

Living Tooth Mould Chart SR Vivodent® S PE

The tooth mould chart contains all the tools necessary for the quick and systematic selection of denture teeth. It includes a shade guide with integrated Facial Meter and Papillameter as well as all the tooth moulds available.

Anterior teeth are arranged according to size and tooth mould. Posterior teeth are assigned to the respective size S, M and L.



Delivery forms

Living Tooth Mould Chart SR Vivodent® S PE	Shade A2	
SR Vivodent S PE / SR Orthotyp S PE	1x	688375