

## BlueLine® Shade Guide

The three layers of original DCL material and the broad range of shades impart a lifelike esthetic appearance to the denture teeth of the BlueLine® Shade Guide.

The shade guide with 16 A–D and 4 Bleach shades is suitable for the shade selection of:

- SR Vivodent® DCL
- SR Vivodent® Plus DCL
- SR Orthotyp® DCL
- SR Ortholingual® DCL
- SR Orthoplane® DCL



### Delivery forms


<b>BlueLine® Shade Guide</b>	<b>1x</b>	<b>613331</b>
------------------------------	-----------	---------------

## SR Vivodent® DCL

The SR Vivodent® anterior tooth line, made of DCL material (Double-Cross-Linked), offers 24 upper and 8 lower anterior A moulds in excellent esthetic quality.

### Advantage:

- Improved abrasion behaviour and easy shade match for combination dentures



### Delivery forms

#### SR Vivodent® DCL Anterior Teeth

16 A–D, 4 Bleach shades  
24 upper, 8 lower

**6 teeth per set**

**564852**

## SR Vivodent® Plus DCL

The SR Vivodent® Plus DCL tooth line comprises the 18 most popular tooth moulds on the North American market.

### Advantage:

- Improved abrasion behaviour and a wide range of upper moulds



### Delivery forms

#### SR Vivodent® Plus DCL Anterior Teeth

16 A–D, 4 Bleach shades  
18 upper

**6 teeth per set**

**591300**