


SR Orthotyp® DCL

Previously well known as the ideal supplement to the SR Vivodent® anterior tooth line, the SR Orthotyp DCL tooth line has been optimized by the use of DCL material (Double-Cross-Linked). The teeth are available in 5 upper and 5 lower posterior tooth moulds.



Advantage:

- Particularly suitable for complete denture as well as for partial denture prosthetics

Delivery forms

SR Orthotyp® DCL Posterior Teeth

16 A–D, 4 Bleach shades		
5 upper, 5 lower	8 teeth per set	565587

SR Ortholingual® DCL

The SR Ortholingual® DCL posterior tooth in 3 upper and 3 lower moulds combines esthetics and functional occlusion with an excellent abrasion behaviour. The SR Ortholingual DCL features an easy and functional setup concept, particularly for lingualized setup.



Advantage:

- Particularly suitable for the lingualized setup method


Delivery forms

SR Ortholingual® DCL Posterior Teeth

16 A–D, 4 Bleach shades		
3 upper, 3 lower	8 teeth per set	565733

SR Orthoplane® DCL

SR Orthoplane® DCL is the 0° posterior tooth available in 3 upper and 3 lower moulds. The most important feature is the simple occlusal anatomy combined with excellent material quality.



Advantage:

- Quick and easy setup
- Particularly suitable for geroprosthodontics

Delivery forms

SR Orthoplane® DCL Posterior Teeth

16 A–D, 4 Bleach shades		
3 upper, 3 lower	8 teeth per set	565843