

Living Tooth Mould Chart BlueLine®

These Living Tooth Mould Charts contain tooth lines in one shade. They are the ideal tool for direct tooth selection – with the patient, but also in the laboratory. All available moulds are arranged according to groups.



Delivery forms

Living Tooth Mould Chart BlueLine®	Shade A2	
SR Vivodent DCL/SR Orthotyp DCL	1x	576324
SR Vivodent DCL/SR Ortholingual DCL	1x	576322
SR Vivodent DCL/SR Orthoplane DCL	1x	576323

SR Vivodent®

The natural appearance and the excellent chemical properties are decisive features of the SR Vivodent® anterior tooth line. The large selection of 20 upper and 8 lower A moulds is based on the tried-and-tested, anatomical Ivoclar Vivadent anterior tooth moulds.



Advantage


- Particularly suitable for complete dentures as well as for partial dentures

Delivery forms

SR Vivodent® Anterior Teeth		
20 Chromascope shades		
20 upper, 8 lower	6 teeth per set	540260

SR Orthotyp®

The 5 upper and 5 lower N moulds of the SR Orthotyp® posterior teeth made of PMMA are the ideal supplement for the SR Vivodent® anterior teeth. They are based on the proven Orthotyp principle.



Advantages:

- Can be combined with SR Vivodent
- The Orthotyp occlusal principle permits quick and efficient tooth setup

Delivery forms

SR Orthotyp® Posterior Teeth		
20 Chromascope shades		
5 upper, 5 lower	8 teeth per set	540050