

Ivostar® | Gnathostar® Shade Guide

The Ivostar® / Gnathostar® shade guide allows the practitioner to determine the tooth shade according to the A–D or Bleach shade system. The shade tabs are made of the original Ivostar® tooth material and are correspondingly layered.



Delivery forms

Ivostar® Gnathostar® Shade Guide	1x	626888
------------------------------------	----	--------

Ivostar®

The attractive and versatile Ivostar® anterior teeth are composed of three layers. They are made of high-quality PMMA, and they are available in esthetic Chromascop, A-D and Bleach shades.



Advantages:

- Large selection of shades
- Harmonized layering
- Versatile application

Delivery forms

Ivostar® Anterior Teeth

20 Chromascop, 16 A-D, 4 Bleach shades 15 upper, 6 lower		
Chromascop shades	sets of 6	517213
A–D and Bleach shades	sets of 6	576694

Gnathostar®

The Gnathostar® posterior teeth are composed of two layers and made of PMMA. The material, layering and shade are optimally matched to the Ivostar® anterior teeth.



Advantages:

- Large selection of shades
- Harmonized layering
- Versatile application

Delivery forms

Gnathostar® Posterior Teeth

20 Chromascop, 16 A-D, 4 Bleach shades 5 upper, 5 lower		
Chromascop shades	sets of 8	517019
A–D shades and Bleach shades	sets of 8	576695