

SR Ivocap® High Impact

The pre-dosed SR Ivocap® High Impact denture base material obtains consistent material quality by the clean, homogeneous mixing of the material in the Cap Vibrator.

The material complies with EN ISO 20795-1.



Advantages:

- Highly impact-resistance and fracture-proof denture base material
- Excellent accuracy of fit without vertical increase
- Outstanding bond with resin teeth
- Optimum polishability
- Tried-and-tested injection technique
- Also polymerizable in the IvoBase Injector

Shades:

Pink; Pink-V; US-L; US-P; S; Preference; US-D

With increased opacity: Pink-V Implant and Preference Implant

Delivery forms

SR Ivocap® Shade Guide	1x	202147
-------------------------------	-----------	---------------

SR Ivocap® High Impact Trial Kit

SR Ivocap High Impact Polymer	10 x 20 g	
Monomer	10 x 30 ml	
Funnels	20x	
Paper inserts	10x	
Pink		531405
Pink-V		531452
Preference		579303
Preference Implant		629096

SR Ivocap® High Impact Standard Kit

SR Ivocap High Impact Polymer	50 x 20 g	
Monomer	50 x 30 ml	
Funnels	100x	
Paper inserts	50x	
Denture bags	100x	
Pink		531396
Pink-V		531451
S		531404
US-L		531402
US-P		531403
Preference		559615
US-D		578936
Pink-V Implant		629097
Preference Implant		629288

Accessories

Plastic floaters, large	1000x	530633
Injection funnels	50x	530615