

ProBase® Cold

The self-curing ProBase® Cold is especially distinguished for its excellent flowing and moulding characteristics. It is suitable for both the pouring as well as the packing technique. Its material properties assure the fabrication of durable denture bases. The material complies with EN ISO 20795-1.



Indications:

- Repairs
- Partial and combination dentures
- Base linings
- Complete dentures

Advantages:

- Convenient processing
- Stability of shape and shade
- For different processing techniques

Shades:

Pink; Pink-V; S; US-L; US-P; 36 P-V; Clear; Preference; US-D
With increased opacity: Pink-V Implant

Delivery forms

ProBase® Shade Guide	1x	531543
ProBase® Cold Trial Kit		
ProBase Cold Polymer	1 x 100 g	
ProBase Cold Monomer	1 x 50 ml	
Dosage system ProBase	1x	
Pink		531486
Pink-V		531487
Clear		531489
ProBase® Cold Standard Kit		
ProBase Cold Polymer	2 x 500 g	
ProBase Cold Monomer	1 x 500 ml	
Dosage system ProBase	1x	
Working accessories:		
Working assortment 1x, mixing cup with cap 1x, mixing spatula 1x		
Pink		531492
Pink-V		531493
US-L		531496
US-P		531497
36 P-V		536281
Clear		531495
Preference		563541
Pink-V Implant		629291
ProBase® Cold Lab Kit		
ProBase Cold Polymer	5 x 500 g	
ProBase Cold Monomer	1 x 1000 ml	
Pink		531498
Pink-V		531499
US-L		531502
US-P		531503
36 P-V		536283
Clear		531501
Preference		578952
US-D		578953

ProBase® Cold Polymer	1 x 500 g	
Pink		531508
Pink-V		531509
S		531510
US-L		531512
US-P		531513
36 P-V		536285
Clear		531511
ProBase® Cold Polymer	2 x 500 g	
Pink		531514
Pink-V		531515
Clear		531517
Pink-V Implant		629294
ProBase® Cold Polymer	5 x 500 g	
Pink		531520
Pink-V		531521
S		531522
36 P-V		536289
Clear		531523
ProBase® Cold Polymer	20 x 500 g	
Pink		531540
Pink-V		531747
US-L		536331
36 P-V		536293
ProBase® Cold Monomer	500 ml	531505
	1000 ml	531506
	4 x 1000 ml	531542
Dosage system ProBase®	1x	531533
Mixing cup / spatula	1x	531534