

SR Triplex® Cold

The self-curing SR Triplex® Cold offers users ease of handling and optimum flow properties.
The material complies with EN ISO 20795-1.



Indications:

- Partial dentures
- Combination dentures
- Base linings
- Repairs

Advantages:

- Convenient processing
- Colours are coordinated with ProBase and SR Ivocap

Shades:

Pink; Pink-V

Delivery forms

SR Triplex® Shade Guide	1x	541894
SR Triplex® Cold Standard Kit		
SR Triplex Cold Polymer	2 x 500 g	
SR Triplex Cold Monomer	1 x 500 ml	
Dosage system SR Triplex	1x	
Pink		541434
Pink-V		541435
36 P-V		547062
SR Triplex® Cold Polymer 2 x 500 g		
Pink		541439
Pink-V		541440
36 P-V		547064
SR Triplex® Cold Polymer 10 x 500 g		
Pink		541441
Pink-V		541442
36 P-V		547065
SR Triplex® Cold Monomer	500 ml	541445
Dosage system SR Triplex®	1x	541654
Plastic container (empty) Polymer for 1000 g	1x	541653