

SR Ivocap® Equipment

The SR Ivocap® system allows the fabrication of fitting methyl methacrylate-based dentures. After mixing, the material is injected into the flask under pressure. Constant pressure (6 bar) during polymerization ensures continuous flow of material. Since the denture base material is continuously injected, the shrinkage of the material is optimally compensated.



Advantages:

- No increase in vertical dimension
- Homogeneous denture bases
- Optimally coordinated components

Delivery forms

SR Ivocap® Equipment		530608
SR Ivocap Pressure Apparatus	1x	
SR Ivocap Clamping Frame	1x	
SR Ivocap Flasks	2x	
SR Ivocap Timer	1x	
SR Ivocap Set of Spring Washers	12x	
Capsule Plunger	1x	
Investment Aids	2x	
SR Ivocap Allen Key	1x	
Accessories		
Pressure apparatus	1x	530610
Clamping frame	1x	530612
Flask	1x	530613
Timer	1x	530614
Set of spring washers	1x	530620

Cap Vibrator

The Cap Vibrator provides reliable and efficient mixing of SR Ivocap® capsules. Its sturdy design and reliable operation ensure homogeneous mixing of pre-dosed materials.



Advantages:

- Pre-dosed material with ideal mixing ratio
- No skin contact during mixing
- Program keys for quick and easy operation
- Convenient capsule holder

Indication:

SR Ivocap capsules

Voltage versions:

220–240 V / 50 Hz, special voltage versions on request

Delivery forms

Cap Vibrator	220–240 V / 50 Hz	530639
Cap Vibrator	1x	
Power cord	1x	