

FH Registration Joint Holder

Used for transferring the UTS 3D registrations to the Stratos® articulators according to the Frankfort Horizontal (FH)



Delivery forms

FH Registration Joint Holder	1x	543468
------------------------------	----	--------

2½ Setting-up Template

For correct tooth set-up of average-valued models



Delivery forms

2½ Setting-up Template	1x	613560
------------------------	----	--------

3D Setting-up Template

For the correct set-up of individually oriented models in conjunction with the UTS® 3D transfer bow, taking the Spee and Wilson curves into consideration.



Important:

This template is not part of the "Personalized" assortment.

Delivery forms

3-D Setting-up Template	1x	536398
-------------------------	----	--------