

## Incisal Plates

L = Incisal plate 15°  
 N = Incisal plate 30°  
 M = Individual incisal plate  
 (Part of the "Personalized" assortment)



### Delivery forms

Incisal Plate		
15°	1x	536391
30°	1x	536392
Individual	5x	536393

## Individual Anterior Guidance Plate

With this accessory, the lateral anterior guidance angle can be individually set to a value between 0° and 55°. The protrusive anterior guidance angle can be infinitely set to a value between 0° and 70°.



### Delivery forms

Individual Anterior Guidance Plate	1x	680584
------------------------------------	----	--------

## Curved Incisal Pin With Holder

This accessory is required for the Individual Anterior Guidance Plate and permits precise, infinite lifting and lowering of the articulator height. A millimetre scale is used as the indicator.



### Delivery forms

Curved Incisal Pin with Holder	1x	680585
--------------------------------	----	--------