

## Model Carriers

- 1) Retentive base block
- 2) Magnetic base block with base former
- 3) Plaster protection plate



### Delivery forms

<b>Magnetic base blocks</b>	2x	536403
	10x	536404
<b>Retentive base blocks</b>	10x	536405
	50x	536406
<b>Plaster protective plates</b>	5x	536408

## Centric Rubber Bands

Centric rubber bands for Stratos 100, 200, 300



### Delivery forms

<b>Centric Rubber Bands</b>	10x	583672
-----------------------------	-----	--------

## Stratos® Facebow Adapter

With this accessory, transfer-/facebows of other manufacturers can be connected with the Stratos articulators (connection pin required).



### Delivery forms

<b>Stratos Facebow Adapter</b>		593851
Left part	1x	
Right part	1x	

## Stratos® Case

Stable case made of aluminium. It enables safe transport of the Stratos articulator and protects it from contamination.



(symbolic figure)

### Advantages

- Stable (aluminium case)
- Compact (fits in standard student lockers)
- Space for Stratos articulator (with model) and accessories
- Including case tag
- Lockable

### Delivery forms

<b>Stratos® Case</b>	1x	583673
----------------------	----	--------