

## SR Nexco®

SR Nexco® is a purely light-curing laboratory composite with micro-opal fillers, suitable for framework-supported and framework-free prosthetic restorations.

In combination with the shades that are tolerant of varying layer thicknesses, true-to-nature optical properties can be achieved for fixed and removable dentures as well as artificial gingiva.



## SR Nexco® Paste Promo Kit



### Delivery forms

#### SR Nexco® Paste Promo Kit A2

645933

SR Nexco Opaquer, shade: A2	1 x 2 ml
SR Nexco Liner, shade: L2	1 x 2 ml
SR Nexco Paste Margin, shade: M2	1 x 2.5 g
SR Nexco Paste Dentin, shade: A2	1 x 2.5 g
SR Nexco Paste Incisal, shade: I2	1 x 2.5 g
SR Gel	1 x 30 ml
SR Link	1 x 5 ml
Various accessories	

#### SR Nexco® Paste Promo Kit A3

667336

SR Nexco Opaquer, shade: A3	1 x 2 ml
SR Nexco Liner, shade: L2	1 x 2 ml
SR Nexco Paste Margin, shade: M2	1 x 2.5 g
SR Nexco Paste Dentin, shade: A3	1 x 2.5 g
SR Nexco Paste Incisal, shade: I3	1 x 2.5 g
SR Gel	1 x 30 ml
SR Link	1 x 5 ml
Various accessories	