

Transil F

Transil F is a light bodied, transparent and therefore light transmitting silicone suitable for use in dental laboratories. Adjusted to the flasking technique, it enables the fabrication of a counter model and therefore the detailed reproduction of the wax-up in the final composite restoration.



Delivery forms

Transil F	2 x 50 ml	673047
incl. Mixing tips	6x	

SR Nexco® Opaquer

The Opaquer paste is distinguished for its high masking power – even in thin coats. Together with SR® Link, it provides for an optimum bond between metal frameworks and the SR Nexco® veneering materials.



Delivery forms

SR Nexco® Opaquer	1 x 2 ml	
A1		640430
A2		640431
A3		640432
A3.5		640433
A4		640434
B1		640435
B2		640436
B3		640437
B4		640438
C1		640439
C2		640440
C3		640441
C4		640442
D2		640443
D3		640444
D4		640445
BL3		640446
BL4		640447

SR Nexco® Liner

The SR Nexco® Liner is used as the basic shade of the prepared tooth in metal-free restorations. It provides a reliable bond between the restoration, luting adhesive and prepared tooth and the SR Nexco® veneering material.



Delivery forms

SR Nexco® Liner	1 x 2 ml	
L1		640476
L2		640477
L3		640478
L4		640479
L5		640480
BL		640481
Clear		640482
Incisal		640483