

SR® Link

SR® Link is a metal-composite bonding agent that provides a covalent bond between the metal framework and resin. It is an easy-to-use and tried-and-tested system.



Delivery forms

SR® Link	1 x 5 ml	545872
----------	----------	--------

SR® Retention Adhesive

This adhesive is a bonding varnish for micro- and macro-retention beads. Its use depends on the space available.



Delivery forms

SR® Retention Adhesive	1 x 20 ml	529498
------------------------	-----------	--------

SR® Micro Retention Beads SR® Macro Retention Beads

Depending on the available space, two sizes of the retention beads are available:

- Micro retention beads 200–300 µm
- Macro retention beads 400–600 µm



Delivery forms

SR® Micro Retention Beads	1 x 15 g	529492
SR® Macro Retention Beads	1 x 15 g	529497

SR® Model Separator

The SR® Model Separator is, among other things, suitable for separating the plaster and the composite veneering materials.



Delivery forms

SR® Model Separator	1 x 10 ml	573105
---------------------	-----------	--------