


OptraSculpt® Pad

OptraSculpt® Pad is a contouring instrument with special foam pad attachments, which is designed for the efficient, non-stick forming and shaping of composites. The highly flexible synthetic foam pads optimally adjust to the anatomical contours and allow a smooth modelling, without leaving any unwanted marks on the composite restoration.



Reference scale on the instrument handle:
The continuous vertical lines show ideal average widths, the dotted vertical lines show average axial inclinations of the maxillary anterior teeth.

Indication:

Suitable for adapting, shaping and contouring uncured sculptable composite materials. It is also suitable for lab composites.

Advantages:

- Non-stick shaping and contouring thanks to foam pad attachments
- Highly elastic foam pads ensure the creation of smooth and even surfaces without unwanted marks
- Reference scales on the instrument handle facilitate professional esthetic results

Delivery forms

| Assortment | | 639791 |
|---|------------|---------------|
| OptraSculpt Pad Instrument | 1x | |
| OptraSculpt Pad Refill/4 mm | 60x | |
| OptraSculpt Pad Refill/6 mm | 60x | |
| OptraSculpt® Pad Refill/4 mm | 60x | 639792 |
| OptraSculpt® Pad Refill/6 mm | 60x | 639793 |
| OptraSculpt® Pad Instrument Refill | 1x | 646004 |

Brush Holder

Disposable brushes and sponges are inserted into the brush holder.



Delivery forms

| | | |
|---------------------------|-----------|---------------|
| Brush Holder white | 1x | 573676 |
|---------------------------|-----------|---------------|