

## Disposable Brushes

The disposable brushes are used to apply liquid materials, such as opaquers, SR® Link, etc.



### Delivery forms

<b>Disposable Brushes</b>	<b>1 x 50</b>	<b>533670</b>
---------------------------	---------------	---------------

## Disposable Sponges

The holders of the disposable sponges were modified in such a way that they fit into the brush holder. The disposable sponges are used to remove the inhibited layer after polymerization.



### Delivery forms

<b>Disposable Sponges</b>	<b>1 x 50</b>	<b>573096</b>
---------------------------	---------------	---------------

## Cannulas

Tips for push syringes (flowable materials) to help achieve more accurate dosing and cleaner application procedures as the extrusion pressure can be controlled.



### Delivery forms

<b>Cannulas</b>	<b>1 x 10</b>	<b>573097</b>
-----------------	---------------	---------------