

## Cannula Caps

These cannula caps prevent the material within the cannula from drying out/polymerizing and protects it from contamination.



### Delivery forms

Cannula Caps	1 x 20	573098
--------------	--------	--------

## SR® Mixing Plate (small)

The SR® Mixing Plate is used to mix the composite components. The light protection prolongs the working time of light-curing materials.



### Delivery forms

SR® Mixing Plate (small)	1x	573100
--------------------------	----	--------

## SR® Mixing Pad (small)

The SR® Mixing Pad is used to avoid having to clean the mixing plate from opaquer residue.



### Delivery forms

SR® Mixing Pad (small)	1x	573099
------------------------	----	--------