

IPS Style® Ceram Intensive Powder Opaquer / IPS Style® Ceram Intensive Paste Opaquer

- **IPS Style® Ceram Opaquer** materials are available in paste and powder form. These opaque ceramic materials have good masking properties for metal frameworks due to their high content of oxyapatite crystals.

Delivery forms

IPS Style® Ceram Paste Opaquer	1 x 5 g
BL1/BL2	701564
BL3/BL4	701565
A1	701566
A2	701567
A3	701568
A3.5	701570
A4	701571
B1	701572
B2	701573
B3	701574
B4	701575
C1	701576
C2	701577
C3	701578
C4	701579
D2	701580
D3	701581
D4	701582
pink	701583

IPS Style® Ceram Intensive Paste Opaquer	1 x 5 g
white	701584
violet	701585
brown	701586
incisal	701587



Delivery forms

IPS Style® Ceram Powder Opaquer 870	1 x 18 g
BL1/BL2	673146
BL3/BL4	673147
A1	673148
A2	673149
A3	673150
A3.5	673151
A4	673152
B1	673153
B2	673154
B3	673155
B4	673156
C1	673157
C2	673158
C3	673159
C4	673160
D2	673161
D3	673162
D4	673163
pink	673164

IPS Style® Ceram Powder Opaquer 870	1 x 80 g
BL1/BL2	673165
BL3/BL4	673166
A1	673167
A2	673168
A3	673169
A3.5	673170
A4	673171
B1	673172
B2	673173
B3	673174
B4	673175
C1	673176
C2	673177
C3	673178
C4	673179
D2	673180
D3	673181
D4	673182
pink	673183



IPS Style® Ceram Intensive Powder Opaquer 870	1 x 18 g
White	673184
Violet	673185
Brown	673186
Incisal	673187