

IPS Style® Ceram Margin

- **IPS Style® Ceram Margin materials** demonstrate a somewhat higher opacity and more fluorescence than the IPS Style Ceram Dentin materials and thus enable the design of ceramic shoulders.
- **IPS Style® Ceram Intensive Margin materials** are used for the individualized design of the ceramic shoulder. They can be mixed with the Margin materials or applied directly.



Delivery forms

IPS Style® Ceram Margin	1 x 20 g
BL	673330
1	673232
2	673233
3	673234
4	673235
5	673236
6	673237
7	673238
IPS Style® Ceram Intensive Margin	1 x 20 g
Yellow	673239
Orange-pink	673240
Orange	673241
Opaque	673242

IPS Style® Ceram Deep Dentin

The IPS Style® Ceram Deep Dentin materials are opaque, shaded dentin materials to be used in areas with limited layer thickness and at incisal ends of frameworks to mask photo-optical “tearing edges” and to achieve true-to-nature results.



Delivery forms

IPS Style® Ceram Deep Dentin	1 x 20 g
BL1 / BL2	673335
BL3 / BL4	673336
A1	673243
A2	673244
A3	673245
A3.5	673246
A4	673247
B1	673248
B2	673249
B3	673250
B4	673251
C1	673252
C2	673253
C3	673254
C4	673255
D2	673256
D3	673257
D4	673258