

IPS Style® Ceram Dentin

The shade and translucency of the IPS Style® Ceram Dentin materials have been adjusted to simulate natural dentin. They permit the correct reproduction of the selected dentin shade.



Delivery forms

IPS Style® Ceram Dentin	1 x 20 g
BL1	673331
BL2	673332
BL3	673333
BL4	673334
A1	673259
A2	673260
A3	673261
A3.5	673262
A4	673263
B1	673264
B2	673265
B3	673266
B4	673267
C1	673268
C2	673269
C3	673270
C4	673271
D2	673272
D3	673273
D4	673274

IPS Style® Ceram Dentin	1 x 100 g
BL1	673339
BL2	673340
BL3	673341
BL4	673342
A1	673275
A2	673276
A3	673277
A3.5	673278
A4	673279
B1	673280
B2	673281
B3	673282
B4	673283
C1	673284
C2	673285
C3	673286
C4	673287
D2	673288
D3	673289
D4	673290

IPS Style® Ceram Mamelon

The IPS Style® Ceram Mamelon materials are intensive, opaquer effect materials for the accentuation of the incisal third. Depending on the working habits of the user, the material is applied in thin stripes on reduced dentin.



Delivery forms

IPS Style® Ceram Mamelon	1 x 20 g
Light	673308
Yellow-orange	673309
Salmon	673310