

IPS Style® Ceram Opal Effect

The IPS Style® Ceram Opal Effect materials are specially shaded incisal materials that permit the simulation of the dynamic photo-optical capabilities of natural teeth.



Delivery forms

IPS Style® Ceram Opal Effect	1 x 20 g
1	673314
2	673315
3	673316
4	673317
5	673318
Violet	673319

IPS Style® Ceram Transpa

The IPS Style® Ceram Transpa materials are available in different shade nuances. They are suitable to reproduce shaded, transparent areas, particularly in the incisal third.



Delivery forms

IPS Style® Ceram Transpa	1 x 20 g
Neutral	673302
Blue	673303
Brown-grey	673304
Orange-grey	673305
IPS Style® Ceram Transpa	1 x 100 g
Neutral	673307

IPS Style® Ceram Special Incisal

The IPS Style® Ceram Special Incisal materials may either be mixed with IPS Style Incisal materials to modify and intensify the shade or applied directly.



Delivery forms

IPS Style® Ceram Special Incisal	1 x 20 g
Yellow	673320
Grey	673321