

## IPS Style® Ceram Incisal Edge

The IPS Style® Ceram Incisal Edge material is used to achieve what is known as the halo effect, which occurs in natural teeth by the refraction of light at the incisal edges.



### Delivery forms

IPS Style® Ceram Incisal Edge	1 x 20 g	673313
-------------------------------	----------	--------

## IPS Style® Ceram Occlusal Dentin

The IPS Style® Ceram Occlusal Dentin materials are used to apply individual characterizations particularly in the occlusal area. They can also be used in the cervical, palatal and lingual areas.



### Delivery forms

IPS Style® Ceram Occlusal Dentin	1 x 20 g	
Orange		673311
Brown		673312

## IPS Style® Ceram Gingiva

The IPS Style® Ceram Gingiva materials are specially shaded ceramic materials for the lifelike reconstruction of gingiva components. They are coordinated with the Gingiva Solutions shade concept from Ivoclar Vivadent and allow true-to-nature esthetics to be achieved due to the comprehensive selection of Gingiva and Intensive Gingiva shades.



### Delivery forms

IPS Style® Ceram Gingiva	1 x 20 g	
1		673344
2		673345
3		673346
4		673347
5		673348
IPS Style® Ceram Intensive Gingiva	1 x 20 g	
1		673349
2		673350
3		673351
4		673352
5		673353
IPS Style® Ceram Basic Gingiva BG34	1 x 20 g	673354