

IPS® Margin Build-Up Liquid

The IPS® Margin Build-Up Liquids are designed for the exclusive use with IPS Style® Ceram. They are available in two variants:

- allround: for producing a consistency that is suitable for conventional processing and drying
- carving: for a shorter processing time and faster drying



Delivery forms

IPS® Margin Build-Up Liquid	1 x 60 ml	
allround		673373
carving		673374

IPS® Build-Up Liquid

The IPS® Build-Up Liquids can be used with IPS e.max® Ceram and IPS Style® Ceram. They are available in two variants:

- allround: for producing a consistency that is suitable for conventional processing and drying
- soft: for a shorter processing time and faster drying



Delivery forms

IPS® Build-Up Liquid	1 x 60 ml	
allround		673375
soft		673376
IPS® Build-Up Liquid	1 x 250 ml	
allround		673377
soft		673378