



IPS d.SIGN®

The fluorapatite leucite glass-ceramic material is distinguished for its lifelike appearance. Fluorapatite is also a component of natural teeth and supports the light optical properties, such as translucency, brightness, and light scattering of the completed restoration.

Advantages:

- Supports the light optical properties
- Opaquer pastes with outstanding masking capabilities

The universal stains and glaze assortment IPS Ivocolor® is used to individually stain and characterize restorations.

Delivery forms

IPS d.SIGN® Opaquer A–D	1 x 3 g	
A1		558198
A2		558199
A3		558200
A3.5		558201
A4		558202
B1		558203
B2		558204
B3		558205
B4		558206
C1		558207
C2		558208
C3		558209
C4		558210
D2		558211
D3		558212
D4		558213
IPS d.SIGN® Intensiv Opaquer	1 x 3 g	
White		556480
Violet		556481
Brown		556482
Incisal		556483