

## IPS d.SIGN® Deep Dentin

With the IPS d.SIGN® Deep Dentin materials, restorations that demonstrate true-to-nature shade effects can be achieved even with limited layer thicknesses. Particularly if space is limited, IPS d.SIGN Deep Dentin materials can be used to intensify the chroma.



### Advantages

- No shade change compared to the basic dentin shade
- Enhanced opacity and shade intensity

### Delivery forms

IPS d.SIGN® Deep Dentin A–D	1 x 20 g	
A1		560512
A2		558258
A3		558259
A3.5		558260
A4		560513
B1		560514
B2		558262
B3		558263
B4		560515
C1		560516
C2		560517
C3		558264
C4		560518
D2 / D3		558261
D4		560519

## IPS d.SIGN® Margin

The IPS d.SIGN® Margin materials permit the design of ceramic shoulders. They are available in Chromascop® and A–D shades. Additionally, there are 4 Intensive Margin materials to achieve special effects. The Margin Add-On material is used for subsequent adjustments of the shoulder area and the corresponding accuracy of fit.



### Advantages:

- High fracture resistance and edge strength
- Excellent firing stability

### Delivery forms

IPS d.SIGN® Margin A–D	1 x 20 g	
A1		560530
A2		558268
A3		558269
A3.5		558270
A4		560531
B1		560532
B2		558272
B3		558273
B4		560533
C1		560534
C2		560535
C3		558274
C4		560536
D2 / D3		558271
D4		560537
IPS d.SIGN® Intensive Margin	1 x 20 g	
Yellow		556562
Orange		556563
Orange-pink		556564
Add-on		556565
Opaque		556566