

IPS d.SIGN® Opaquer Liquid

This liquid is used to slightly dilute dried opaquer. It must not be used for mixing opaquers in paste form.



Delivery forms

IPS d.SIGN® Opaquer Liquid	1 x 15 ml	556645
----------------------------	-----------	--------

IPS d.SIGN® Margin Build-Up Liquid

Used to mix the IPS d.SIGN® Margin materials



Delivery forms

IPS d.SIGN® Margin Build-Up Liquid	1 x 60 ml	556648
------------------------------------	-----------	--------

IPS d.SIGN® Build-Up Liquid Medium

This liquid is excellently suitable for mixing IPS d.SIGN® materials for most layering techniques and restorations. The IPS d.SIGN Build-Up Liquid Medium is particularly suitable for technicians who prefer a longer working time and enhanced stability of IPS d.SIGN. In this way, IPS d.SIGN may be processed in smaller portions, without the constant need of moistening. Given the longer working time than with the Build-Up Liquid Light, a longer closing time (pre-drying time) is required for the firing procedure.

Delivery forms

IPS d.SIGN® Build-Up Liquid Medium	1 x 60 ml	556644
	1 x 250 ml	558985

