

IPS d.SIGN® Build-Up Liquid *Light*

This liquid is excellently suitable for mixing IPS d.SIGN® materials used for smaller restorations or for smaller portions in conjunction with the 2nd dentin and incisal firing (corrective firing). The Build-Up Liquid Light is suitable for technicians who prefer drier working procedures, i.e. without the constant need to suction off moisture, since the liquid evaporates more quickly. In this way, IPS d.SIGN can be conveniently applied with an instrument.



Delivery forms

IPS d.SIGN® Build-Up Liquid Light

1 x 60 ml	564045
1 x 250 ml	564043

IPS d.SIGN® Build-Up Liquid *Premium*

This liquid is excellently suitable for mixing IPS d.SIGN® material used for larger restorations, for which longer plasticity and smooth consistency is desired. The Build-Up Liquid Premium is suitable for technicians who prefer a long working time and moist material coupled with high stability. In this way, IPS d.SIGN can be processed without constant moistening. Since the material stays moist for a prolonged period of time, however, controlled suction is required, together with a longer closing time (pre-drying time) for the firing procedure.



Delivery forms

IPS d.SIGN® Build-Up Liquid Premium

1 x 60 ml	564046
1 x 250 ml	564044