

IPS d.SIGN® Build-Up Liquid Optimum 1

Excellently suitable for mixing IPS d.SIGN®, with its outstanding firing behaviour and lowest possible shrinkage. The materials are easily suctioned off and condensed and are particularly suitable for the manual reduction of materials with the help of instruments.



Delivery forms

IPS d.SIGN® Build-Up Liquid Optimum 1

	1 x 60 ml	575541
	1 x 250 ml	575542

IPS d.SIGN® Build-Up Liquid Optimum 2

Excellently suitable for mixing IPS d.SIGN®, with its properties that support the build-up, as well as the optimized firing properties with low shrinkage. This liquid is particularly suitable for targeted material application without continuous suction.



Delivery forms

IPS d.SIGN® Build-Up Liquid Optimum 2

	1 x 60 ml	575543
	1 x 250 ml	575545