

IPS InLine® PoM Ingots

With only seven shades, all Chromascop®, A–D and the new BL Bleach shades can be easily reproduced. Consequently, it is possible to press over different patient cases in one press cycle. Thanks to two ingot sizes, there isn't unnecessary material waste. The final tooth shade is achieved by individual characterization with the IPS InLine® Shade/Stains and Glaze materials.



Delivery forms

IPS InLine® PoM Ingots S	1 x 5
BL	602394
1	602395
2	602396
3	602397
4	602398
5	602399
6	602400

IPS InLine® PoM Touch Up

The seven IPS InLine® Touch Up materials, which are especially coordinated with the ingot shades, are ideal for final shape and shade adjustments.



Delivery forms

IPS InLine® PoM Touch Up	1 x 20 g
BL	602401
1	602402
2	602403
3	602404
4	602405
5	602406
6	602407