

## IPS InLine® System Build-Up Liquid L

The liquid is suitable to mix IPS InLine® ceramic materials for smaller restorations or for smaller portions. The Build-Up Liquid L is particularly suitable for dental technicians who prefer a drier working consistency – i.e. without continuous suction – since the liquid evaporates much quicker. In this way, IPS InLine can be easily manipulated with an instrument.



### Delivery forms

IPS InLine® System Build-Up Liquid L	1 x 60 ml	593348
	1 x 250 ml	593351

## IPS InLine® System Build-Up Liquid P

The liquid is suitable for mixing IPS InLine® ceramic materials that will be used for fabricating larger restorations where a longer working time and smooth consistency is desired. The Build-Up Liquid P is suitable for technicians who prefer a long working time and moist material coupled with high stability.



### Delivery forms

IPS InLine® System Build-Up Liquid P	1 x 60 ml	593349
	1 x 250 ml	593352