

## IPS® Model Sealer

The IPS® Model Sealer is used to seal the plaster model. The IPS Model Sealer closes the pores of the plaster and prevents the moisture in the ceramic material from being absorbed by the plaster.



### Delivery forms

IPS® Model Sealer	1 x 50 ml	593401
-------------------	-----------	--------

## IPS® Margin Sealer

The IPS® Margin Sealer Liquid is used to seal the plaster die when fabricating a ceramic shoulder.



### Delivery forms

IPS® Margin Sealer	1 x 20 ml	593402
--------------------	-----------	--------

## IPS® Ceramic Separating Liquid

This Ceramic Separating Liquid is used on areas that have been sealed with either the IPS® Model Sealer or the IPS® Margin Sealer Liquid. Use of this liquid will ensure a clean separation between the ceramic material and the plaster model and/or die.



### Delivery forms

IPS® Ceramic Separating Liquid	1 x 15 ml	593403
--------------------------------	-----------	--------