

IPS® Shade V

The ready-mixed IPS® Shade dentin stains in paste form permit subsequent modification of shades on the fired ceramic restoration. Additionally, they can be mixed with the ceramic materials for individual characterization.



Advantages:

- Easy, time-saving handling
- High resistance to temperature changes and oral conditions

Delivery forms

IPS® Shade V	1 x 3 g
A1	531777
A2 / A3 / A3.5	531778
A4	531779
B1	531780
B2 / B3 / B4	531781
C1 / C2	531782
C3 / C4	531783
D2 / D3	531784
D4	531785

IPS® Stains-P

These stains in paste form permit easy, individualized characterization of the restorations. The stains may be either directly applied on the ceramic surface or mixed with the ceramic materials. The countless possibilities of mixing the shades ensure virtually unlimited creativity.



Advantages:

- High resistance to temperature changes and oral conditions
- Broad temperature range

Delivery forms

IPS® Stains-P	1 x 3 g
White	531590
Bamboo-beige	531591
Caramel-brown	531592
Copper-brown	531593
Kork-brown	531594
Mahogany	531595
Azure	531596
Black	531597
Orange	531598
Basic yellow	531599
Basic red	531600
Basic blue	531601