

## IPS Classic® Correction Powder

Low-fusing add-on material for IPS Classic® and IPS Classic V metal-ceramics. Particularly suitable for small corrections and shape adjustments of fired restorations.



### Advantages:

- Low firing temperature of 690°C / 1274°F
- Mixing with IPS Shade or IPS Stains-P possible
- High resistance to acid and oral conditions

### Delivery forms

IPS Classic® Correction Powder	1 x 20 g	541671
--------------------------------	----------	--------

## IPS Classic® Glaze Paste

Easy-to-use glazing material in paste form for a true-to-nature gloss on ceramic restorations.



### Delivery forms

IPS Classic® Glaze Paste	1 x 3 g	531602
--------------------------	---------	--------