

## IPS Classic® Build-Up Liquid L (Long-term)

Build-Up Liquid "L" has been developed for technicians who require a longer working time, as well as for markets in warmer regions. The liquid contains an additive that raises the boiling point and ensures slower evaporation of the liquid. This results in a prolonged working time while the excellent build-up properties remain unchanged.



### Delivery forms

IPS Classic® Build-Up Liquid L	1 x 60 ml	529398
	1 x 250 ml	529402

## IPS Classic® Build-Up Liquid N (Standard)

This build-up liquid is particularly suitable for technicians who apply the ceramic materials in small portions using a brush. Build-Up Liquid "N" provides outstanding stability and demonstrates optimum build-up properties. This material is used by technicians who prefer moist conditions, i.e. no excessive suction. Moreover, Build-Up Liquid "N" is suitable for everyday use and for most of the layering techniques.



### Delivery forms

IPS Classic® Build-Up Liquid N	1 x 60 ml	529397
	1 x 250 ml	529401

## IPS Classic® Build-Up Liquid S (Spatula Technique)

This liquid is particularly suitable for technicians who frequently use a suction device, condense the materials and subsequently reduce the build-up using an instrument. The build-up liquid is easily sucked off the contoured crown. In this way, a relatively solid mass is achieved that can be individually shaped using an instrument.



### Delivery forms

IPS Classic® Build-Up Liquid S	1 x 60 ml	529399
	1 x 250 ml	529403