

## IPS e.max® Press Multi (Polychromatic) Ingots

The IPS e.max® Press Multi ingots are available in Bleach BL and A–D shades. IPS e.max Press Multi ingots are used for pressing highly esthetic veneers, anterior and posterior crowns and hybrid abutment crowns quickly and efficiently. The restorations show a shade and translucency progression similar to that of natural teeth. The ingots are suitable to be processed in the staining technique.



### Delivery forms

IPS e.max® Press Multi	1 x 5 each	
BL2		638638
A1		638036
A2		638147
A3		638148
A3.5		638149
B1		638151
B2		638152
C1		638155
C2		638156
D2		638159

## IPS e.max® Press HT (High Translucency) Ingots

Given their high translucency, the HT ingots are ideally suitable for being pressed to full contour (staining technique). They are available in 16 A–D and 4 Bleach BL shades.



### Delivery forms

IPS e.max® Press HT	1 x 5 each	
BL1		626352
BL2		626353
BL3		626354
BL4		626355
A1		626320
A2		626321
A3		626322
A3.5		626323
A4		626324
B1		626325
B2		626326
B3		626327
B4		626328
C1		626329
C2		626330
C3		626331
C4		626332
D2		626333
D3		626334
D4		626335

IPS e.max® Press HT L	1 x 3 each	
BL1		626356
BL2		626357
BL3		626358
BL4		626359
A1		626336
A2		626337
A3		626338
A3.5		626339
A4		626340
B1		626341
B2		626342
B3		626343
B4		626344
C1		626345
C2		626346
C3		626347
C4		626348
D2		626349
D3		626350
D4		626351