

Multilink® Hybrid Abutment

Multilink® Abutment is a self-curing luting composite for the permanent cementation of lithium disilicate (LS₂) glass-ceramic and zirconium oxide ceramic structures on implant sleeves made of titanium/titanium alloy and zirconium oxide.



Indication:

Extraoral cementation of lithium disilicate (LS₂) glass-ceramic and zirconium oxide ceramic structures (abutment and abutment crowns) on implant sleeves made of titanium/titanium alloy and zirconium oxide.

Advantages:

- High bonding strength for a reliable cementation
- Easy handling with the practical Automix syringe

Shades:

HO 0

Delivery forms

Multilink® Hybrid Abutment Starter Kit

662810

Multilink Hybrid Abutment Syringe HO 0	1 x 2.5 g
Mixing Tips	4x
Monobond Plus	1 x 1 g
Mixing Pad	1x
Various accessories	

Refill

Multilink Hybrid Abutment Syringe	1 x 9 g
Mixing Tips	15x
Various accessories	
HO 0	647182

IPS e.max® Abutment Solutions Cem Kit

638959

Multilink Hybrid Abutment Syringe HO 0	1 x 9 g
Mixing Tips	15x
Virtual Extra Light Body fast	1 x 50 ml
Mixing Tips	48x
Intra Oral Tips	100x
Monobond Plus	1 x 5 g
Liquid Strip	1 x 2.5 g
Various accessories	