

IPS e.max® CAD Crystallization Tray

The IPS e.max® CAD Crystallization Tray is a firing tray for the crystallization of IPS e.max CAD restorations as well as for the Stains and/or Glaze firing of IPS e.max ZirCAD and IPS Empress CAD restorations with the IPS e.max CAD Crystall./materials. To conduct the firing, the restorations are placed on the IPS e.max CAD Crystallization Pins supplied with the tray.



Delivery forms

IPS e.max® CAD Crystallization Tray	605367	
IPS e.max CAD Crystallisation Tray	1x	
IPS e.max CAD Crystallization Pin (sizes: S, M, L)	9x	

IPS e.max® CAD Speed Crystallization Tray

IPS e.max® CAD Speed Crystallization Tray is a special firing tray designed for the speed crystallization of IPS e.max CAD restorations and for the Stains and/or Glaze firing of IPS e.max ZirCAD restorations with IPS e.max Crystall./materials.




Delivery forms

IPS e.max® CAD Speed Crystallization Tray	647633	
IPS e.max CAD Speed Crystallization Tray	1x	
IPS e.max CAD Speed Crystallization Pin	3x	

IPS e.max® CAD Crystallization Pins

The pins assist you with the positioning of restorations on the IPS e.max® CAD Crystallization Tray or the IPS e.max CAD Speed Crystallization Tray.



Delivery forms

IPS e.max® CAD Crystallization Pin Refill		
S / M / L	6 x each	605368
XS	3 x	656105