

IPS e.max® Ceram Impulse Occlusal Dentin

The Occlusal Dentin materials are available in two different shades. They are used as the shade basis for the occlusal surface and to intensify the chroma.



Delivery forms

IPS e.max® Ceram Impulse Occlusal Dentin	1 x 20 g each	
Orange		596983
Brown		596984

IPS e.max® Ceram Impulse Mamelon

The Mamelon materials are available in three different shade gradations. They demonstrate high opacity and excellent masking capabilities even when applied in very thin layers. Depending on the working habits of the user, these materials are applied in thin lines on the reduced dentin. In this way, an individualized appearance of the incisal third is achieved.



Delivery forms

IPS e.max® Ceram Impulse Mamelon	1 x 20 g each	
Light		596986
Yellow-orange		596987
Salmon		596988

IPS e.max® Ceram Impulse Incisal Edge

This material is used to achieve what is known as the “halo effect”, which is caused in natural teeth by the light refraction at the incisal edge.



Delivery forms

IPS e.max® Ceram Impulse Incisal Edge	1 x 20 g	596996